



NLM Malaria Screener User Guideline

Outline

- Steps for Collaboration Page 3
- Malaria Screener User Manual Page 4
- Frequently Asked Questions Page 17

Steps for Collaboration




1. Get equipment

- Microscope with 100X objective lens
- Android smartphone with 12+ Mega Pixel camera.
- Adapter

Recommendation: <https://www.telescopeadapters.com/128-universal-mounts>


* Refer to Question 1 in Frequently Asked Questions for more details on equipment



2. Email us

- Email us a Gmail address to join the beta test group

hang.yu@nih.gov



3. Take images

- Uses the app to take images of 10+ thin or thick smears with manual counts.




*No personal data will be collected

Malaria Screener User Manual




Proper Smear Imaging

To achieve optimal results, proper images should be presented to the app.




- Center Position, field-of-view should fit to image



- Even illumination, avoid shadow (mostly along with border)




- Clean lens




Proper Smear Imaging


- Cells distribution




Well separated single cells if possible



Too dense





Too sparse




Too sparse




- Sharp image, no blur



Blood film diagnosis




NLM Malaria Screener



Start a New Session


Database

- ① On the main page of the app, press "Start a New Session".



- ② Select smear type.

Blood film diagnosis



Here, thin smear is chosen as an example to show the workflow of a session.

Blood film diagnosis

Search patient ID to fill out the form automatically if the patient is in the database

Fill out the forms of patient and slide information in the next 3 pages.

6 Input patient information.

Fields with asterisk marks are required fields, other fields can be left empty. (Only Patient ID and Slide ID are required fields)

7 Fields with asterisk marks are required fields, other fields can be left empty. (Only Patient ID and Slide ID are required fields)

Adjust white balance when necessary to produce image with true to life color

Smaller value gives a classifier that is less sensitive to infected cells, and vice versa.


You can also change the amount of cells you want to collect during each session.

Settings

Database



- ① On the main page of the app, press “Database”.




- ② Data for all the patients are saved in the database.




- ③ Slides for each patient

Database




- ④ Click “VIEW IMAGES” to view all the images of the slide.



Tap the image to view it in full screen

- ⑤ Each section displays the original and result images side.




- ⑥ View image in full screen


Recording manual counts

It is very important for us to get the manual counts of each image so that we can evaluate the performance of our app. You can either input the manual counts during an imaging session or afterwards. Please continue to read the next few slides for the detailed instructions.

Enter manual counts during a session




- 1 Find the more options icon on the top right corner of the result page and press "Enter Manual Counts".



- 2 Enter your manual counts in the dialog that appears in the center of the page

Enter manual counts in the database


Ideally, manual counts should be entered during the screening session as shown in the previous page. This is only an alternative method.



- 1 Find the more options icon on the top right corner of the "Slide Image" page and press "Enter Manual Counts".



- 2 The manual counts for each image are listed on this page. Press each list item to edit.



- 3 Another option is to long press on the image thumbnail, and enter your counts in the dialog that appears.

Recording manual counts


The app can export its database to an csv file. Information of all the saved sessions are recorded in this file including the manual counts that the user have entered. When you upload your images, this file will be automatically included in the folder.

	A	B	C	D	E	F	G
1	patient_id	slide_id	image_id	cell_count	infected_count	cell_count_gt	infected_count_gt
2	17	NLMM17	20180702_152506.png	234	5	N/A	13
3	17	NLMM17	20180702_152659.png	224	7	N/A	13


Manual count of total number of RBCs. (Can be omitted if it's too much work.)

Manual count of Infected RBCs.


Upload data



- ① On the main page of the app, press “database”.



- ② Press “UPLOAD”.




- ③ Press “YES” and enter registration page.




***User needs a Box account in order to use this feature. You can register for a free account at [Box.com](https://www.box.com).**

Upload data



A screenshot of a mobile application's 'Register' screen. The screen has a blue header with a back arrow and the word 'Register'. Below the header are several text input fields: 'Last Name *', 'Email Address *', 'Organization *', 'Department *', 'Address line 1 *', 'Address line 2', 'City *', 'Zip Code *', and 'Country *'. At the bottom of the form is a blue button labeled 'REGISTER'.

4 Fill out the form and then press “REGISTER”.





A screenshot of a mobile application's login screen for Box. The screen has a blue header with the 'box' logo. Below the header is a large blue icon representing a folder. The main text reads 'Log in to grant access to Box'. There are two input fields: 'Email Address' and 'Password'. Below these fields is a blue button labeled 'Authorize'. Underneath the button are links for 'Use Single Sign On (SSO)' and 'Forgot password'. At the bottom, there is a small line of text: 'By granting NLM_MalariaScreener_Upload_UserAuth access to Box, you are agreeing to Box's Terms of Service and Privacy Policy.' The footer shows '©2020 Box'.



5 User will be prompted to enter credentials to grant access to the account.

Upload data



*The app will only perform uploading tasks with a WiFi connection to avoid potential data charges.



- ⑥ **One-click upload.** Simply press the “Quick Upload” button to upload all images that have not been uploaded.





- ⑦ **Upload with manual selection.** User can also manually select the folders.



- ⑧ Upon uploading, a floating widget indicates the upload progress.




NLM Malaria Screener




- ⑨ **Single Session Upload.** After each session, the app will attempt to upload data from the session. This is to split up the workload.

Image acquisition mode



- 1 At the beginning of a session, turn on Image Acquisition Mode in Settings.



- 2 Find ideal field and capture images like usual. The image will not be processed, and the app will immediately go to the next page when an image is accepted.



- 3 The counts will all be zeros since no processing is done in this mode. Press "End Session" when you decide enough images have been captured.

Frequently Asked Questions

1. What smartphone and/or microscope adapter should I order?

Answer: Many smartphones and microscope adapters are good to be used for our app. The key here is for the two to work well with each other, however this is difficult to tell without testing. We give a few tips that will help you when choosing your devices:

A few tips while choosing smartphones & microscope adapters:

- Processor power of the smartphone will determine the processing speed, so a newer model is usually preferred.
- The buttons on the sides of smartphone shall not clash with the clamps of the adapters.
- Smartphones with multi-lens camera are not recommended.
- Go with good quality adapters. A cheap microscope adapter may work. However, it may cause you a lot of trouble.

For example,

- a. It maybe be difficult to get it attached to the microscope.
- b. It maybe be difficult to align the camera with the eyepiece.

Recommendations:

Smartphone: Samsung Galaxy S Series + Microscope adapter: <https://www.telescopeadapters.com/128-universal-mounts>

Frequently Asked Questions

2. Does the app collect patient's personal information?

Answer: No, the app does not collect any patient's personal information.



U.S. National Library of Medicine



National Institute of
Allergy and
Infectious Diseases



U.S. National Library of Medicine
Lister Hill National Center for Biomedical Communications